

Sonometer (Brass)



Product Categories: [Physics Equipments](#), [Sound And Waves Oscillations](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/sonometer-brass/>

Product Description

A Sonometer brass is a device for demonstrating the relationship between the frequency of the sound produced by a plucked string, and the tension, length and mass per unit length of the string. It is Mostly used in School, college, research, university Laboratories.

Feature:

Sonometer brass Soft wood nicely polished with metal fitting.

soft wood With wooden riders of sonometer

And provide with 3 different wire 100cm long.

Brass fitting without wt.

Brass fitting with Teak wood box.