

## Sonometer (Softwood)



**Product Categories:** [Physics Equipments](#), [Sound And Waves Oscillations](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/sonometer-softwood/>

### Product Description

A Sonometer brass is a device for demonstrating the relationship between the frequency of the sound produced by a plucked string, and the tension, length and mass per unit length of the string. It is most commonly used in all kinds of Laboratories.

**Available in size** - 6"x 6" (size box)

### Feature

Sonometer Softwood is made of nicely polished.

With brass fitting 6"x 6" box size.

Tuning fork set of 8

Tuning fork set cheap quality

Accessories for experiment

Metal set of 8 Welch type

Raggs eng land (Imported set of 8)

Rubber pad

Spare knife edge