

Spare Boss Head



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/spare-boss-head/>

Product Description

Spare Boss Head - Laboratory Instrument

Boss Head- Made of cast iron hold upto 13 mm dia T type screw powder coating Made of cast iron hold upto 18 mm dia T type screw powder coating. Made of cast iron hold upto 25 mm dia T type screw powder coating.

Features:-

BOSS HEAD: Made of cast iron hold upto 13mm dia T Type screw powder coating.

BOSS HEAD : Made of Al-alloy die casting mild steel screw upto 13mm dia powder coated finish.

Oxidized with boss head.

Moulded round (stable painted).

Moulded C.P