

## Spare Die & Punches



**Product Categories:** [Pathology Lab](#)

**Product Page:** <https://www.labappara.com/product/spare-die-punches/>

### Product Description

#### **Spare Die & Punches - Laboratory Equipment**

Spare die & punches of chromium steel.

Available in size - 6, 9, & 12 mm (any one size).

Uniformity across all punches and dies giving interchangeability between every punch or die within the same set.

Punch height measurement using highly accurate programmable height gauge.

Computerized heat treatment.

Optimum clearance prevents product loss from lower punches.

#### **Features:**

Robust construction

Accurate dimensions

High strength

Anti-corrosive

Precisely designed

Flawless finish