

Spare Perimeter Tubes



Product Categories: [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/spare-perimeter-tubes/>

Product Description

Spare Polarimeter Tubes - Physics Optical Equipment

Extra corning glass brass collars with windows, bubble trap in the long. These tubes are used in Polarimeter.

Rotating Triple Shadow effect with fine optics, suitable for sugar industries without Lamp Assembly

Polarimeter Tubes 100 mm / 200 mm / 400 mm - Center / Side Filling