


Spectrometer Prism (Dense Flint)



Product Categories: [Educational equipments](#), [Glass Prism](#), [Light & Optics](#), [Optical Instruments](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/spectrometer-prism-dense-flint/>

Product Description

Spectrometer Prism (Dense Flint) - Physics Optical Equipment

Made from Borosilicate glass, these equilateral prisms are optically worked with 2 faces polished and most suitable for use with a spectrometer.

Optically worked for use with Spectrometer, Equilateral, two faces polished 32 x 32 mm height. Optically worked with refractive index of 1.62

Equilateral

2 sides polished

Dent flint glass

Size:-

(a) 25x 25mm (Refractive index)-1.62

(b) 32x 32mm

(c) 38x 38mm