


## Spectrometer Prism (Extra Dense Flint)



**Product Categories:** [Educational equipments](#), [Glass Prism](#), [Light & Optics](#), [Optical Instruments](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:**

<https://www.labappara.com/product/spectrometer-prism-extra-dense-flint/>

### Product Description

Spectrometer Prism (Extra Dense Flint) - Physics optical Equipment

Made from Borosilicate glass, these equilateral prisms are optically worked with 2 faces polished and most suitable for use with a spectrometer.

Optically worked with refractive index of 1.65

Equilateral,

2 sides polished

Extra dent flint glass

**Size:-** (a) 25x 25mm      (Refractive index)-1.65

(b) 32x 32mm

(c) 38x 38mm