

# SPECTROPHOTOMETER



**Product Categories:** [Analytical Lab Instruments](#), [Spectrophotometer](#)  
**Product Tags:** [Analytical Laboratory Equipment](#)  
**Product Page:** <https://www.labappara.com/product/spectrophotometer/>

## Product Description

SPECTROPHOTOMETER - Analytical Laboratory Equipment

A spectrometer is an instrument that measures the amount of photons (the intensity of light) absorbed after it passes through sample solution.

**Features:**

Rugged structure  
Resistance to corrosion  
Accurate reading  
Water proof

**Details:**

Cover 195 to 1000nm  
4 Position Cuvette Holder  
RS-232C Computer Interface  
High wavelength Precision  
Highly Accurate & Stable