

## Spirit Lamp Stainless Steel (40z)



**Product Categories:** [Chemistry Equipments](#), [General Chemistry Labware](#)

**Product Tags:** [chemistry lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/spirit-lamp-stainless-steel-40z/>

### Product Description

Spirit Lamp Stainless Steel (40z)- Scientific Equipment

A Spirit Lamp is laboratory equipment used to produce an open flame.

It is widely used in chemistry laboratory to perform experiments

Spirit lamp available with us are of brass, glass, stainless steel or aluminium.

Stainless Steel Spirit Lamp Sheet die pressed, with woven wick in metal holder, screw on cap.

Capacity 100ml.