

Spirit Level



Product Categories: [Chemistry Equipments](#), [General Chemistry Labware](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/spirit-level/>

Product Description

Spirit Level - Scientific Equipment

A spirit level, bubble level or simply a level is an instrument designed to indicate whether a surface is horizontal (level) or vertical (plumb).

Available in size - 2", 3", 4", 6"

AVAILABILITY:

With plastic.

With all brass.

With aluminum.