

Spray Tower



Product Categories: [Engineering Equipment](#), [Mass Transfer Lab](#)

Product Page: <https://www.labappara.com/product/spray-tower/>

Product Description

Spray Tower

Spray Tower setup is designed to demonstrate principles of spray extraction and mass transfer operation in a column.

The set up consists of a glass column, feed tanks, air regulator with pressure gauge, receiving trays and piping. Rota-meter are provided for flow measurement

Features:

Stainless Steel tanks and wetted parts

Special arrangement for changing height of interface zone

Continuous operation

Superb Painted structure

Simple to operate & maintain

Compact & stand alone set up