

SPRING BALANCE



Product Categories: [Balance & Scale](#), [Chemistry Equipments](#), [Measurements](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/spring-balance-2/>

Product Description

Spring Balance - Physics Instrument

Flat form, high impact polystyrene body zero adjustment having anodized metal scale with clear and accurate division. High tension special wire spring is used for more durability with ring & hook packed in box. Color coded spring balances for both force & mass measurement. Individually calibrated dual scale spring Balances has zero adjustment device with stainless steel suspension loop and load hook.

Capacity - 100gm, 250gm, 500gm, 1000gm