

SPRING BALANCE (Transparent)



Product Categories: [Balance & Scale](#), [Chemistry Equipments](#), [Measurements](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/spring-balance-transparent/>

Product Description

Spring Balance - Physics Instrument

Transparent Acrylic with ring for suspension and hook for heavy weight. The Spring Balance has a tubular metal body with suspension ring and load hook and with clear scaling. Moreover, the Heavy Duty Spring Balance are available in the array and provide accurate result in respective applications.

Capacity - 100gm, 250gm, 500gm, 1000gm, 2kg, 5kg.