

## Stage Micrometer



**Product Categories:** [Microscope](#), [Microscope Assecories](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/stage-micrometer/>

### Product Description

#### Stage Micrometer - Laboratory Equipment

A Stage Micrometer is simply a microscope slide with a finely divided scale marked on the surface. The scale is of a known true length and is used for calibration of optical systems with eyepiece articulate patterns.

#### **Features:**

Precision Glass Stage Micrometer Slide

Least Count 10 Microns

Scale of 1 mm divided in 100 Parts

High Contrast Sharp graduated lines

Fine ground Glass

Ideal for Microscope, camera & software Calibration

Comes in protective case