

Stalganometer



Product Categories: [Laboratory Glassware](#), [Viscometer Borosilicate Glass](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/stalganometer/>

Product Description

Stalganometer - Laboratory Equipment

The stalagmometric method is one of the most common methods for measuring surface tension. A stalagmometer is a device for investigating surface tension using the stalagmometric method. It is also called a stactometer or stalogometer. Straight, graduated.

Features:

Precisely designed

Perfect readings

Easy to clean