

## Stream Boiler (Copper)



**Product Categories:** [Heat And Thermodynamics](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/steam-boiler-copper/>

### Product Description

Steam Boiler - Physics Equipment

A boiler is a closed vessel in which water or other fluid is heated.

Heavy thick body made of copper.

Use to generate steam

Heat energy is used to convert water into steam

This steam may be used for different purposes.

Including water heating, central heating, boiler-based power generation, cooking, and sanitation.

5 LITER

1 LITER

2 LITER