

Stereo Binocular Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/stereo-binocular-microscope/>

Product Description

Stereo Binocular Microscope - Biology Laboratory Equipment

Available in different configurations for specific AP/publications, with dual illumination, for surface and transmitted illumination, versatile stand can be used for specimen up to 60 mm thickness.

Focusing: Vertical movement controlled by coarse focusing large sized knobs. Rack and Pinion focusing system, single stroke up to 50 mm.

Eyepieces: 10 x (paired) with extra wide field.

Objectives: Paired Achromatic Objectives 1 x & 3 x or 2 x & 4 x.

Illumination: Built-in illumination for incident (EPEE) light and sub-stage (Transmitted) light.

Features:

Durable finish standards

Precise design

Temperature resistance

Easy to use

Lightweight