

Stomach Model



Product Categories: [Anatomical Models](#), [Biology Lab](#), [Educational equipments](#), [Educational Models](#), [Organ Models](#)

Product Tags: [biology lab](#), [fibre model](#), [science model](#)

Product Page: <https://www.labappara.com/product/stomach-model/>

Product Description

Stomach Model- Educational Model

Stomach Model- Stomach Model With the longitudinal section, the model shows the structures of the gastric folds, pyloric valve, pyloric sphincter muscles, gastric mucous and the transitional mucous of the gastric-esophagus. Made of plastic and magnified 2 times of the natural size.