

## Storage Boxes



**Product Categories:** [Laboratory Plastic Ware](#), [Storage Boxes Carriers](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/storage-boxes/>

### Product Description

#### Storage Boxes - Scientific Laboratory Equipment

These storage boxes, made of polypropylene, prove to be an excellent option for storage of different powdered reagents & other material in the laboratory. Square shape of these boxes help in saving shelf space. Available in three different sizes, these boxes have polypropylene caps, are autoclavable, provide good chemical resistance & heavy duty service.

#### **Available in sizes :**

250 ml

500 ml

1000 ml