

## Student Microscope



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:** <https://www.labappara.com/product/student-microscope/>

### Product Description

#### Student Microscope - Laboratory Equipment

Commonly used where high magnification is not required for like in school labs . heavy die cast body with white lens bright field fixed condenser, two eyepiece 10 x, 15 x, two objectives 10 x, 45 x (SL). Well packed in thick Plywood box.

#### **Features:**

Clear and sharp view

Compact size

Easy to use

Precisely engineered

Long functional life

Require less maintenance

#### **Specifications:**

Student Polarising Microscope: Having circular stage roatable in 360 degree,

graduate with

verier reading upto 0.1

Objectives: 5 x, 10 x and 45 x

Eye Pieces: WF 10 x and H 15 x