

Tablet Coating Pan



Product Categories: [Pathology Lab](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/tablet-coating-fan/>

Product Description

Tablet Coating Pan - Laboratory Equipment

12" dia., stainless steel, the machine is worked by oil immersed gear box with F.H.P. motor drive. Air blower is operated independently and air is thermostatically controlled to give temperature 30 to 100 degree C. Switches for motor, blower, heater and thermostat are fixed on a panel with neon indicators.

Specifications :

MS Sheet, Powder coated body, Stainless Steel Pan

Provided with main lead 15 amp plug

Operated by geared motor that works through timer and air blower, which is further operated by energy regulator

Working of the machine through oil immersed gearbox with $\frac{1}{4}$ H.P Motor resulting or by heavy duty DC geared motor

Tablet coating pan 12" dia stainless steel

Air blower is operated independently and air is thermostatically controlled.

Tablet coating Pan to give Temperature: 30o C TO 110o C

Switches for blower & motor are fixed on a panel through neon indicators