

Tablet Counting Device



Product Categories: [Pathology Lab](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/tablet-counting-device/>

Product Description

Tablet Counting Device - Laboratory Equipment

Tablet counting device for 100 tablets of 9.5mm.

Perforations made of aluminum sheet.

Tablet counting device provided with pouring lip and handle.

CAPACITY :

100 tablets

250 tablets

500 tablets

Capsules & Oval Tablets, 100