

## Tablet Hardness Tester (Monsanto type)



**Product Categories:** [Pathology Lab](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:**

<https://www.labappara.com/product/tablet-hardness-tester-monsanto-type/>

### Product Description

Tablet Hardness Tester (Monsanto type) - Laboratory Equipment

Tablet hardness tester (Monsanto type) is a laboratory technician used by the laboratory to test the breaking point and structural integrity of a tablet under conditions of storage, and handling before usage.

The breaking point of a tablet is based on its shape.

Monsanto type small portable device to measure the hardness.

#### **Monsanto type:**

A fairly instrument for taking quick reading for the hardness of Tablets

Made using brass turned parts

Duly chromium plated with scale length of 66 mm

Calibrated to read between 0 to 20 kg of hardness division

The unit is enclosed in plastic case

#### **Application:**

Used to determine the edge-to-edge (diametric) tablet hardness

## **Features:**

Superior functionality

High accuracy

Easy to operate