

Tablet Hardness Tester (Monsanto type)



Product Categories: [Pathology Lab](#)

Product Tags: [Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/tablet-hardness-tester-monsanto-type/>

Product Description

Tablet Hardness Tester (Monsanto type) - Laboratory Equipment

Tablet hardness tester (Monsanto type) is a laboratory technician used by the laboratory to test the breaking point and structural integrity of a tablet under conditions of storage, and handling before usage.

The breaking point of a tablet is based on its shape.

Monsanto type small portable device to measure the hardness.

Monsanto type:

A fairly instrument for taking quick reading for the hardness of Tablets

Made using brass turned parts

Duly chromium plated with scale length of 66 mm

Calibrated to read between 0 to 20 kg of hardness division

The unit is enclosed in plastic case

Application:

Used to determine the edge-to-edge (diametric) tablet hardness

Features:

Superior functionality

High accuracy

Easy to operate