

## Tablet Hardness Tester (Peizer type)



**Product Categories:** [Pathology Lab](#)

**Product Tags:** [Scientific Equipment](#)

**Product Page:**

<https://www.labappara.com/product/tablet-hardness-tester-peizer-type/>

### Product Description

Tablet Hardness Tester (Peizer type) - Scientific Equipment

Tablet hardness testing, is a laboratory technique used by the laboratory to test the breaking point and structural integrity of a tablet under conditions of storage, and handling before usage.

The breaking point of a tablet is based on its shape.

#### **Pfizer type:**

Suitable for tablet hardness testing

Very handy and portable

Records hardness quickly and directly on dial indicator

When a tablet is crushed by pressing the jaws of the instrument, the force is recorded on a dial indicator calibrated in force unit 1/5 kg

Squeeze of hand grip generates the braking force

Release of the hand pressure opens the jaws

Accuracy :  $\pm 0.2$ kg Least Count : 0.2Kg

Safe testing Force : 20kg

Maximum Tablet Diameter : 25 mm

**Features:**

Superior functionality

High accuracy

Easy to operate

**Capacity - 20 kg**