

## Tablet Making Machine



**Product Categories:** [Pharmaceutical Equipments](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/tablet-making-machine/>

### Product Description

Tablet Making Machine - Pharmaceutical Equipment

Automatic compression and ejection device with each revolution of the wheel drive.

Automatic feed and removal of tablets.

Heavy cast body with assentric pressure supreme.

Stainless steel hopper provided and connected by a rubber tube.

Easy setting for thickness and hardness of tablets with nuts & locknut and set of two spanners.

Die and set punches are of chromium steel fitted on stand.

### Feature

Tablet Making Machine with consists of hand operated.

Tablet Making Machine with consists of electrically operated.

Single punch for making tablets up to 1000 mg as per size.

Tablet Making Machine of the dia and punch are available are also 6, 9, 12 mm.

Easily rotating wheel for making tablets

Compressing lever for adjustable harness of tablets

Adjustable lower punch for adjusting the thickness

Stainless steel day and punch size 6, 9, 12 mm.

Electrical model can be provided in same with ½ h.p motor.