


Teclu Burner



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/teclu-burner/>

Product Description

Teclu Burner - Scientific Equipment

The Teclu burner is a laboratory gas burner, a variant of the bunsen burner. It can produce a hotter flame than a bunsen burner. Teclu Burners : Nickel Plated Brass Burner Tube with rotatable Air Regulator and cylindrical rifflled connector, Burner produces a very hot flame. Suitable for glass work and soldering. elegantly powder coated Aluminum Base.

Feature -

Duco painted.

On cast iron base.

Wide & stable base provide.

Packed in box.

Burner tube 60 × 10 mm (height × diameter)

Base 80 mm diameter

Overall Height 95 mm

Inlet tube 10 mm mean OD