

## Telescopic Delivery Table



**Product Categories:** [Hospital Equipment & Furniture](#)

**Product Tags:** [hospital](#), [medical](#)

**Product Page:** <https://www.labappara.com/product/telescopic-delivery-table/>

### Product Description

Telescopic Delivery Table - Hospital Equipment

Table top in three section backrest adjustable on ratchet. Middle section with U cut & leg section can slide completely under main frame. The whole top can be trendelenburg position with the help of gas spring.

The table completely with following accessories .

Railing on head and side, hand, grip, saline stand, stainless steel bow and lithotomic with strip.

**Overall Size** - 1870 L X 680 W X 800 H