

Test Tube Baskets



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/test-tube-baskets/>

Product Description

Test Tube Baskets - Scientific Laboratory Equipment

Used For storing, cleaning, and autoclaving test tubes and other laboratory glassware These autoclavable, unbreakable and non-corrosive Test Tube Basket are ideal for rinsing and draining glassware. Tapered sides for easy handling and compact storage. The basket's cage like construction ensures the safe storage of test tubes or culture tubes during sterilization or incubation.

Particulars	Packing
-------------	---------

140 x 120 x 110 mm	6 Pcs
--------------------	-------

160 x 160 x 160 mm	6 Pcs
--------------------	-------