

Test Tube Brush



Product Categories: [Chemistry Equipments](#), [General Chemistry Labware](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/test-tube-brush/>

Product Description

Test Tube Brush - Scientific Equipment

A test tube or spout brush is a brush used for cleaning test tubes and narrow mouth laboratory glassware, such as beakers and flasks.

It is composed of nylon, synthetic, or animal fur bristles of various diameters lined against a rather sturdy wire handle with a looped end for hanging.

The function of a test tube brush is to clean test tubes used in a science laboratory.

Shape Fan-shaped end