

Test Tube Stand (Round)



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [Laboratory Equipments](#), [Laboratory plastic ware](#)

Product Page: <https://www.labappara.com/product/test-tube-stand-round/>

Product Description

Test Tube Stand (Round) - Scientific Laboratory Equipment

A handy, little, space saver, this circular twelve place Polypropylene Test Tube Stand has four holes & eight holes in its top plate. The base has matching hemispherical wells & vertical pins for drying tubes surrounded by a trough to catch drippings. This autoclavable rack can be easily disassembled for cleaning purpose. This auto-clavablerack can be easily disassembled for cleaning purpose.

Material: Polypropylene

Available in: 19 mm & 25 mm Tubes