


Test Tube Stand (Wire Pattern)



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/test-tube-stand-wire-pattern/>

Product Description

Test Tube Stand (Wire Pattern) - Scientific Laboratory Equipment

An economical substitute for wire racks, here are these submersible, autoclavable racks which are stackable when empty. These are available in different colors and are made of a special blend of Polypropylene that makes it sink in a water bath and maintain their stability even in agitated water. These space saving racks can easily be disassembled for cleaning purpose. the base is in the form of a grid with square openings and fold up sides and end bars that lock into rack tops.

Particulars	Packing
13 mm X 36 Tubes	4 Pcs
16 mm X 36 Tubes	4 Pcs