


Test Tube Stand



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/test-tube-stand/>

Product Description

Test Tube Stand - Scientific Laboratory Plastic ware Equipment

This autoclavable Stand is perhaps the most commonly used Test Tube Stand in various laboratories.

For better visibility, this place rack is end supported by columns.

The top plate has holes for holding test tubes, four holes to hold tubes where as the base has corresponding hemispherical well.

It also has six vertical pins fro drying tubes.

Another option available in this section is a test tube stand in which all the holes in the top plates can hold tubes.

The top plate has 2 holes for holding 25 mm dia test tubes, four holes to hold 16 mm dia tubes where as the base has corresponding hemispherical wells.

Particulars	Packing
16 mm & 25 mm	12 Pcs
25 mm	12 Pcs