


Test Tube Stands:- 6 Hole Test Tube Stand



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [Lab Measuring Instruments](#), [Laboratory Equipments](#)

Product Page:

<https://www.labappara.com/product/test-tube-stands-6-hole-test-tube-stand/>

Product Description

Test Tube Stands:- 6 Hole Test Tube Stand - Scientific Laboratory Plastic ware
Test tube stand is a laboratory equipment that is used to hold upright multiple test tubes at the same time. It is most commonly used when various different solutions are needed to work with simultaneously, for safety reasons, for safe storage of test tubes, and to ease the transport of multiple tubes. The top plate has 2 holes for holding 25 mm dia test tubes, four holes to hold 16 mm dia tubes where as the base has corresponding hemispherical wells. It also has six vertical pins for drying tubes.

Particulars	Packing
16 mm & 25 mm	12 Pcs
25 mm	12 Pcs