

## Test Tubes



**Product Categories:** [Laboratory Glassware](#), [Test Tubes / Culture Tubes](#), [Test Tubes / Culture Tubes Borosilicate Glass](#)

**Product Tags:** [glassware](#), [medical](#)

**Product Page:** <https://www.labappara.com/product/test-tubes-2/>

### Product Description

Test Tubes - Laboratory Glassware Equipment

"JSGW", with Inter changeable Stopper

#### Features

made of borosilicate glass 3.3

very good chemical resistance

minimal thermal expansion

high resistance to thermal shocks

with white marking area

with screw thread and black screw cap including Butyl/PTFE sealing insert

sterilisable up to 110 °C

wall thickness: approx. 1.0 to 1.1 mm

ungraduated

Available in capacity - 5 ml, 10 ml, 20 ml, 25 ml, 50 ml, 100 ml, 150 ml.