

Test Tubes Round



Product Categories: [Laboratory Glassware](#), [Test Tubes / Culture Tubes](#), [Test Tubes / Culture Tubes Borosilicate Glass](#)

Product Tags: [glassware](#)

Product Page: <https://www.labappara.com/product/test-tubes-round/>

Product Description

Test Tubes Round - Laboratory Glassware Equipment

Made of glass medium walled with rim for use in school, colleges, universities, research industrial laboratories.

Wall thickness 1 mm comes in pack of 100 pieces.

Test Tubes Round made from Borosilicate 3.3 exp. glass tubing, Round bottom, Rimless, Plain.

Available in capacity - 3 ml, 5 ml, 7 ml, 10 ml, 15 ml, 25 ml, 35 ml, 50 ml, 75 ml, 100 ml, 175 ml