

Tetrode/ Pentode Valve Characteristics Apparatus



Product Categories: [E/M Apparatus & Valve Experiments](#), [E/M Apparatus & Valve Experiments](#), [Electronics](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/tetrode-pentode-valve-characteristics-apparatus/>

Product Description

Tetrode/ Pentode Valve Characteristics Apparatus

Features :

Instrument comprises of Two DC Regulated Power Supplies 0-300VDC/30mA & One 0-15 VDC/10mA, AC Power Supply 6.3 VAC/1Amp, four round meters for voltage & Current measurement, valve diagram is Printed on Front panel, connections of Supplies & Valve brought out at 4mm Socket/Terminals.