


# Time Division Multiplexer Trainer



**Product Categories:** [Basic Communication Lab](#), [Engineering Equipment](#)

**Product Page:** <https://www.labappara.com/product/time-division-multiplexer-trainer/>

## Product Description

Time Division Multiplexer Trainer  
Built in +/-5V/350mA DC Power Supply  
Data Generator 1 K b/s/Channel  
Low High, Logic Low sources  
Time Division Multiplexer - 8 Ch  
Time Division Demultiplexer - 8 Ch  
LED indication for channel outputs  
Assembled in Wooden Box with circuit  
screen printed on Glass Epoxy PCB with  
2 mm socket for test points & to see the waveforms  
Set of Patch Chords & Experimental  
Manual