

Tincture Press



Product Categories: [Pathology Lab](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/tincture-press/>

Product Description

Tincture Press- Laboratory Equipment

Extra heavy for making tinctures decoctions & infusions.

Tincture Press Stainless steel disc with screw drive.

Tincture Press a stainless steel inner perforated vessel.

Following are the features of these machines:

Additional heavy duty for making tinctures

Decotious infusion cap

Disc made up of pure stainless steel

Drive for fixing screws

Mesh made from stainless steel in the form of an inner drilled container with stop cock

Capacity – 1 liter