

## Triode Valve Characteristics Apparatus



**Product Categories:** [E/M Apparatus & Valve Experiments](#), [E/M Apparatus & Valve Experiments](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:**

<https://www.labappara.com/product/triode-valve-characteristics-apparatus/>

### Product Description

#### Triode Valve Characteristics Apparatus

**Features :**

Instrument comprises of DC Regulated Power Supply 0-300VDC/ 30mA & 0-25VDC/10mA, AC Power Supply 6.3 VAC/1Amp, two round meters for voltages & current measurement, valve diagram is printed on front panel, connections of supplies & valve brought out at 4mm socket.