

Tripod Stand (Wire type)



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/tripod-stand-wire-type/>

Product Description

Tripod Stand (Wire type) -

The tripods are made of Mild steel (no Cast iron).

Tripods are with Nylon feet.

Joints are nicely argon welded.

Feet are covered with rubber stumps.

A tripod is a portable three-legged frame, used as a platform for supporting the weight and maintaining the stability of some other object.

A tripod provides stability against downward forces and horizontal forces and movements about horizontal axes.

A laboratory tripod is a piece of three-legged equipment commonly used to conduct experiments in laboratories.

Available in size-

Wire type & Moulded regular quality

A

6"X 3"

B

8"X 3"