

Tubular Steel Trolley



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [biology lab equipment](#), [Hospital Equipment](#)

Product Page: <https://www.labappara.com/product/tubular-steel-trolley/>

Product Description

Tubular Steel Trolley - Hospital Equipment

Strong tubular steel frame with six smooth running drawers.

Mobile on 15cm swivel castors.

Provide with one SS IV rod , massage board, six handout bins & cylinder cage.

Fixed with buffer rail.

Pre treated, epoxy coated finishing.

Overall Size - 100W X 85D X 152H cm