

Tuning Forks



Product Categories: [Physics Equipments](#), [Sound And Waves Oscillations](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/tuning-forks/>

Product Description

Tuning Forks - Physics Equipment

Tuning fork use for producing sound waves by vibrating .

The fork consists of a handle and two vibrating metal strip.

When tuning fork hit with rubber hammer, the two vibrating metal strip begin to vibrate.

The vibration of tuning folk produce disturbances of surrounding air molecules.

Feature -

Tuning folk set consist of different frequency Set of 8.

Frequencies of Tuning folk starts from 256 to 512 steel C.P.

High tensile strength

Accuracy in pouch

Excellent performance