

Tunnel Diode Characteristics Apparatus

Product Categories: [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

Product Page:

<https://www.labappara.com/product/tunnel-diode-characteristics-apparatus/>

Product Description

Tunnel Diode Characteristics Apparatus

Technical Specifications :

Inbuilt Fixed DC regulated power supply

Output Voltage : +5VDC

On Board Digital Panel Meter

Digital panel meter for measuring

V1 : Voltage across Resistance R 3

V2 : Voltage across Tunnel Diode

Potentiometer & Diode

Potentiometer : R1 (Current control)

Tunnel Diode : IN 3717

High quality bakelite used as front panel & mounted on light weight shock proof plastic cabinet

Symbol diagram printed on bakelite Front Panel & all important test Points are brought out on front panel

Power requirement : 230 VAC 10%, 50Hz.

Weight : 1.5Kg Approx.

Standard Accessories : Power chord, Patch chords & Instruction manual.