

Tunning Forks (Aluminum Packed)



Product Categories: [Physics Equipments](#), [Sound And Waves Oscillations](#)

Product Tags: [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/tunning-forks-aluminum-packed/>

Product Description

Tunning Fork (Aluminium Packed) - Physics Instument

Tuning fork aluminum use for producing sound waves by vibrating.

The fork consists of a handle and two vibrating aluminum strip.

When tuning fork hit with rubber hammer, the two vibrating metal strip begin to vibrate.

The vibration of tuning folk produce disturbances of surrounding air molecules

Forks are individually constructed of Aluminum.

Feature -

Tuning fork of Made of aluminum packed in wooden box

Set of 8 in a wooden box.

High tensile strength

Accuracy in pouch

Excellent performance