

Tunning Forks (Blue Steel)



Product Categories: [Physics Equipments](#), [Sound And Waves Oscillations](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/tunning-forks-blue-steel/>

Product Description

Tunning Forks - Physics Instrument

Tuning fork use for producing sound waves by vibrating.

The fork consists of a handle and two vibrating strip. When tuning fork hit with rubber hammer, the two vibrating metal strip begin to vibrate.

The vibration of tuning folk produce disturbances of surrounding air molecules.

Forks are individually constructed of blue steel.

Feature -

Tuning fork of Made of blue steel packed in wooden box

Set of 8 in a blue steel .