

Tunning Forks (Chromium plated)



Product Categories: [Physics Equipments](#), [Sound And Waves Oscillations](#)

Product Tags: [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/tunning-forks-chromium-plated/>

Product Description

Tunning Forks (Chromium plated) - Physics Equipment

Tuning fork Chromium plated use for producing sound waves by vibrating .

The fork consists of a handle and two vibrating metal strip. When tuning fork hit with rubber hammer, the two vibrating metal strip begin to vibrate.

The vibration of tuning folk produce disturbances of surrounding air molecules.

Feature -

The turning folk is consists of Pye type.

The turning folk is Chromium plated for long life.

Set of 8 in velvet case.

Bright chrome plated.