

Tunning Forks (Welch type)



Product Categories: [Physics Equipments](#), [Sound And Waves Oscillations](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/tunning-forks-welch-type/>

Product Description

Tuning fork welch type use for producing sound waves by vibrating.

The fork consists of a handle and two vibrating Welch strip. When tuning fork hit with rubber hammer, the two vibrating metal strip begin to vibrate.

The vibration of tuning folk produce disturbances of surrounding air molecules

Forks are individually constructed of tuning fork in Welch type.

Feature -

Tuning fork of Made of welch type packed in velvet box.

Set of 8 in a velvet box.

Turning forks is large size in chrome plated of packed in velvet box.

And then welch type special quality tested very accurate in velvet box.

It Is used in School, College, University, Research and Other Laboratories.