

## Tunning Forks (Wooden)



**Product Categories:** [Physics Equipments](#), [Sound And Waves Oscillations](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/tunning-forks-wooden/>

### Product Description

#### **Tunning Forks (Wooden) - Physics Instruments**

Tuning fork Wooden use for producing sound waves by vibrating .

The fork consists of a handle and two vibrating wooden strip. When tuning fork hit with rubber hammer, the two vibrating wooden strip begin to vibrate.

The vibration of tuning folk produce disturbances of surrounding air molecules

#### **Feature -**

Tuning fork set of 8 in a wooden case.

Tuning forks are individually constructed of forged wooden.

Accurate & Mild steel C.P.

Tuning forks are individually constructed of forget steel.

Tuning fork set of 8 in a wooden case.