

## UJT Characteristics Apparatus



**Product Categories:** [Characteristics And Application](#), [Electronics](#)

**Product Page:** <https://www.labappara.com/product/ujt-characteristics-apparatus/>

### Product Description

UJT Characteristics Apparatus **Objective :** To plot V-I Characteristics & to measure inter-base resistance.

**Features :** Instrument comprises of Two DC Regulated Power supplies 0-15VDC/150mA & 0-30VDC/150mA, three round meters for voltage & current measurement, one UJT 2N2646 mounted behind the panel, connections of Supplies, Meters & UJT brought out at 4mm Sockets.